



WebDisplay Install Guide

WebDisplay Installation/ Upgrade

- 1 From your Internet browser, log onto the MedPACS Web page: www.medpacs.com
- 2 At the "Home" page click-on "Support"
- 3 At the "Support" page, under "Upgrades" click-on "WebDisplay Upgrade"
- 4 At the "File Download" screen click on the "Open" button (this will download the WebDisplay viewer onto your system – depending on the quality of your connection the file transfer may take some time)
- 5 At the "Welcome" screen (WebDisplay Setup program) click on "Next"
- 6 At the "Setup Complete" screen, click-on "Finish"
- 7 Plug the WebDisplay Key into the USB port on the computer. If the LED light on the Key doesn't light, reboot your computer (you can leave the Key plugged-in). If after rebooting the LED is still not on, contact MedPACS.
- 8 Contact MedPACS* for WinRSF setup and installation test. Note that you will not be able to receive images until this setup and installation test is completed by MedPACS; after a successful test, you'll be ready to view images.



** To minimize support charges, call during normal business hours.*

*The above set-up procedure needs to be completed at each computer/
workstation that the WebDisplay will be used at to view images*



Ph: 262-367-0181
Fx: 262-367-9574
support@medpacs.com